

DIGIPYRO PROJECT

Installing DigiPyRo

Samuel May

June 30, 2016

1 Mac (OSX 10.10 or later)

1. Install **Xcode & the Command Line Tools** from the App Store (this may take 10-20 minutes). Then open **Xcode**, agree to the terms and open a terminal to install the command line tools. Enter the following in your terminal:

```
xcode-select --install
```

2. Install **Homebrew**. This is a package manager, which allows you to easily install different packages that we need from the command line. Open a Terminal and enter the following:

```
cd ~
```

```
ruby -e "$(curl -fsSL https://raw.githubusercontent.com/Homebrew/install/master/install)"
```

For additional information on **Xcode** and **Homebrew**, see:

<https://coolestguidesontheplanet.com/installing-homebrew-on-os-x-el-capitan-10-11-package-manager-for-unix-apps/>

3. Install **Python2** (with the help of **Homebrew**!). In your Terminal enter the following:

```
brew update
```

```
brew install python
```

Note: It is important that you install **Python2** (which is installed by the command `brew install python`)! **DigiPyRo** will not work with **Python3**, so you must follow this step even if you have already installed **Python3**.

4. Install **OpenCV**. Enter the following in your terminal:

```
brew tap homebrew/science
```

```
brew install opencv
```

installing **OpenCV** may take 10-20 minutes.

5. Link **Python** and **OpenCV** together. Enter the following in your terminal:

```
cd ~
```

```
touch .bash_profile (this creates the file if it does not already exist)
```

```
cat ~/.bash_profile | grep PYTHONPATH
ln -s /usr/local/Cellar/opencv/2.4.13/lib/python2.7/site-packages/cv.py cv.py
ln -s /usr/local/Cellar/opencv/2.4.13/lib/python2.7/site-packages/cv2.so cv2.so
```

Note: the most recent releases of Python and OpenCV may have changed since this was written. You may need to adjust your version numbers accordingly:

```
ln -s /usr/local/Cellar/opencv/2.?.?.lib/python2.?.site-packages/cv.py cv.py
ln -s /usr/local/Cellar/opencv/2.?.?.lib/python2.?.site-packages/cv2.so cv2.so
```

In order to find your version numbers, simply navigate to the appropriate directory and type the command:

```
ls
```

to display the contents of the directory. You can then read your version number and adjust the previous two symlink commands accordingly.

Note: the `ln -s` command creates what is called a “symlink”. Basically, your computer will not know where to search for Python and OpenCV, so we facilitate this process by creating aliases in your home directory.

6. Install the necessary Python modules. We will use `pip`, Python’s built-in package manager which stands for the recursive title “`pip` installs packages”. Enter the following in a terminal:

```
pip install numpy
pip install matplotlib
pip install scipy
```

7. Download DigiPyRo. Enter in a terminal:

```
git clone https://github.com/sam-may/DigiPyRo/tree/master
```

8. Check if everything is working! First let’s check if Python and OpenCV are correctly linked. Open Python by entering the following in a terminal:

```
python
```

then enter:

```
import cv2
```

in the Python interpreter. If you see no errors then Python and OpenCV are correctly linked!


Now, let’s try running DigiPyRo. Open a terminal and enter:

```
cd ~
```

```
cd DigiPyRo
```

```
python DigiPyRo.py
```

If this menu appears:



The screenshot shows a window titled "DigiPyRo" with the following controls:

- Curvature of table (in RPM, enter 0 for a flat surface): 0.0
- Physical rotation (of camera, RPM): 0.0
- Additional digital rotation (RPM): 0.0
- ☐ Custom-Shaped Mask (checking this box allows for a polygon-shaped mask. default is circular)
- Full filepath to movie: [empty text box]
- Save output video as: [empty text box]
- Start and end times (in seconds): 0.0 [empty text box] 0.0
- Frames per second of video: 0.0
- ☐ Add distance units calibration
- Length unit (e.g. cm, ft): [empty text box]
- Unit count: 0.0
- ☐ Track Ball
- ☐ Create plots with tracking results
- [Start!] button

then you have successfully installed DigiPyRo! Well done.

2 Windows

http://opencv-python-tutroals.readthedocs.io/en/latest/py_tutorials/py_setup/py_setup_in_windows/py_setup_in_windows.html